

SAY GOODBYE TO FINGER PRICKING AND RETIRE YOUR TEST STRIPS TOO.



“ NO NEEDLE, NO BLOOD, NO TEST STRIPS, NO PAIN, NO INFECTION, HUGE SAVINGS”



 **Bluetooth®**



GET BLOOD GLUCOSE READING ANYTIME , ANYWHERE WITHOUT ANY PAIN.

WORLD'S FIRST PATENTED INDUCTIVE NON-INVASIVE BLOOD GLUCOSE METER.

TAIWAN , USA , UK , JAPAN , GERMANY , EUIPO , CHINA PATENT RIGHT , KOREA , INDIA PATENT PENDING.

ETM-G01E Non-invasive Blood Glucose Meter

Features:

- * Compliant to test 3 body modes for diabetes: Normal, Pre-Diabetes & Type 2 Diabetes.
- * Put thumbs pulp and forefinger pulp of the right and left hands on the 4 sensors within 10 seconds reaction time.
- * Easy operation (just turn on the power, select user mode then touch the test points).
- * Easy maintain (no need to change the needle, no need to clean the test parts).
- * Easy carry (no clean kits, no blood collection pen, no alcohol or cotton required).
- * Memory capacity 1,000 sets test result.
- * One product for your whole family.
- * Power source: 2 * "AAA" batteries.
- * Bluetooth built-in.
- * Size : 96 * 86 * 24 (mm).
- * Designed and made in Taiwan.

mg/dL (mmol/L)	FBG Fasting Blood Glucose	PBG Postprandial Blood Glucose
Normal	80 ~ 100 (4.4 ~ 5.6)	120 ~ 140 (6.7 ~ 7.8)
Pre-Diabetes	101 ~ 125 (5.6 ~ 6.9)	141 ~ 200 (7.8 ~ 11.1)
Type 2 Diabetes	126 + (7.0 +)	201 + (11.2 +)